

PRP-X200 Analytical Guard Cartridge Starter Kit (1 holder, 2 cartridges), Stainless Steel

Specification Sheet (Part/REF # 79456)



PRP-X200 is a versatile type of cation exchange resin particularly adept at interacting with both inorganic and organic cations, including alkali and alkaline earth metals. It is capable of handling a variety of cations, whether they are monovalent or divalent, providing wide-ranging applications in several scientific and industrial processes.

Product Specifications

Separation Mode	Cation Exchange
Stationary Phase	PRP-X200
Format	Guard Cartridge Kit
Hardware Inner Diameter	2.1 mm
Hardware Length	20 mm
Particle Size	10 μm
Packing Material Type	PSDVB/Sulfonic Acid
Particle Surface Area (m ² /cm ³)	450
Pore Volume (cm ³ /g)	0.752
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Organic/Inorganic Cations
UOM	1 holder, 2 cartridges

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79456>

Spec. Version B