

PRP-X800 Analytical Guard Column Starter Kit (1 holder, 2 cartridges), Stainless Steel

Specification Sheet (Part/REF # 79830)



Protect your PRP-X800 column.

Product Specifications

Separation Mode	Cation Exchange
Stationary Phase	PRP-X800
Format	Guard Cartridge Kit
Hardware Inner Diameter	2.1 mm
Hardware Length	20 mm
Particle Size	7 μm
Packing Material Type	PSDVB/Itaconate
Particle Surface Area (m ² /cm ³)	50
Pore Volume (cm ³ /g)	0.166
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Mono-Divalent Organic/Inorganic Cations - Transition Metals
UOM	1 holder, 2 cartridges

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79830>

Spec. Version B