

PRP-X800 Analytical Guard Column Replacement Cartridges (5/pk), PEEK

Specification Sheet (Part/REF # 79833)



Protect your PRP-X800 column.

Product Specifications

Separation Mode	Cation Exchange
Stationary Phase	PRP-X800
Format	Replacement Guard Cartridge
Hardware Inner Diameter	3.0 mm
Hardware Length	8 mm
Particle Size	7 µm
Packing Material Type	PSDVB/Itaconate
Guard Column Size	Analytical
Hardware Material	PEEK
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Mono-Divalent Organic/Inorganic Cations - Transition Metals
Frit Material	Titanium
UOM	5/pk

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79833>

Spec. Version B