

PRP-X300 Analytical Guard Cartridge Starter Kit (1 holder, 2 cartridges), PEEK

Specification Sheet (Part/REF # 79373)



Protect your PRP-X300 column.

Product Specifications

Separation Mode	Ion Exclusion
Stationary Phase	PRP-X300
Format	Guard Cartridge Kit
Hardware Inner Diameter	3.0 mm
Hardware Length	8 mm
Particle Size	12-20 µm
Packing Material Type	PSDVB/Sulfonic Acid
Guard Column Size	Analytical
Hardware Material	PEEK
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Organic Acids/Alcohols
Frit Material	Titanium
UOM	1 holder, 2 cartridges

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79373>