

# PRP-X300 Analytical Guard Cartridge Replacement Cartridges (5/pk), Stainless Steel

Specification Sheet (Part/REF # 79453)



Protect your PRP-X300 column.

## Product Specifications

<b>Separation Mode</b>	Ion Exclusion
<b>Stationary Phase</b>	PRP-X300
<b>Format</b>	Replacement Guard Cartridge
<b>Hardware Inner Diameter</b>	2.1 mm
<b>Hardware Length</b>	20 mm
<b>Particle Size</b>	12-20 $\mu$ m
<b>Packing Material Type</b>	PSDVB/Sulfonic Acid
<b>Guard Column Size</b>	Analytical
<b>Hardware Material</b>	Stainless Steel
<b>Max Pressure</b>	5,000 psi
<b>Mobile Phase Limits</b>	pH 1-13. 0-100% aqueous, organic modifier
<b>Compound Classification</b>	Organic Acids/Alcohols
<b>UOM</b>	5/pk

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79453>