

PRP-X300 Semiprep/Preparative Guard Cartridge Starter Kit (1 holder, 1 cartridge), Stainless Steel

Specification Sheet (Part/REF # 79129)



Protect your PRP-X300 column.

Product Specifications

Separation Mode	Ion Exclusion
Stationary Phase	PRP-X300
Format	Guard Cartridge Kit
Hardware Inner Diameter	4.6 mm
Hardware Length	20 mm
Particle Size	12-20 µm
Packing Material Type	PSDVB/Sulfonic Acid
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Exchange Capacity	0.17 meq/gm
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Temp Limits	5-60°C
Compound Classification	Organic Acids/Alcohols
UOM	1 holder, 1 cartridge

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79129>

