

PRP-h5 Analytical Guard Column Starter Kit (1 holder, 2 cartridges), Stainless Steel

Specification Sheet (Part/REF # 79267)



Protect your PRP-h5 column.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-h5
Format	Guard Cartridge Kit
Hardware Inner Diameter	2.1 mm
Hardware Length	20 mm
Particle Size	10 μm
Packing Material Type	PSDVB
Particle Surface Area (m^2/cm^3)	277
Pore Volume (cm^3/g)	0.526
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Peptides/Proteins
UOM	1 holder, 2 cartridges

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79267>

Spec. Version B