

PRP-X400 Analytical Guard Column Starter Kit (1 holder, 2 cartridges), Stainless Steel

Specification Sheet (Part/REF # 79224)



Protect your PRP-X400 column.

Product Specifications

Separation Mode	Cation Exchange
Stationary Phase	PRP-X400
Format	Guard Cartridge Kit
Hardware Inner Diameter	2.1 mm
Hardware Length	20 mm
Particle Size	7 µm
Packing Material Type	PSDVB/Sulfonic Acid
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Organic/Inorganic Cations
UOM	1 holder, 2 cartridges

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79224>

Spec. Version B