

PRP-X400 Analytical Replacement Cartridges (5/pk), Stainless Steel

Specification Sheet (Part/REF # 79225)



Protect your PRP-X400 column.

Product Specifications

Separation Mode	Cation Exchange
Stationary Phase	PRP-X400
Format	Replacement Guard Cartridge
Hardware Inner Diameter	2.0 mm
Hardware Length	20 mm
Particle Size	7 μ m
Packing Material Type	PSDVB/Sulfonic Acid
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Compound Classification	Organic/Inorganic Cations
UOM	5/pk

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79225>

Spec. Version B