

PRP-X500 Analytical Guard Column Replacement Cartridges (5/pk), PEEK

Specification Sheet (Part/REF # 79320)



Protect your PRP-X500 column.

Product Specifications

Separation Mode	Anion Exchange
Stationary Phase	PRP-X500
Format	Replacement Guard Cartridge
Hardware Inner Diameter	3.0 mm
Hardware Length	8 mm
Particle Size	12-20 µm
Packing Material Type	Poly(methacryl amido-propyl trimethyl ammonium chloride)
Particle Surface Area (m ² /cm ³)	1.4
Pore Volume (cm ³ /g)	0.0063
Guard Column Size	Analytical
Hardware Material	PEEK
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Large Proteins/Labeled DNA
Frit Material	Titanium
UOM	5/pk

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79320>

Spec. Version B