

PRP-1 Semiprep/Preparative Guard Column Starter Kit (1 holder, 1 cartridge), Stainless Steel

Specification Sheet (Part/REF # 79121)



Organic compounds: small molecules (< 2,000 mw), oligonucleotides, pharmaceuticals, steroids, nucleic acids, vitamins, herbicides.

Product Specifications

| | |
|--------------------------------|---|
| Separation Mode | Reversed-Phase |
| Stationary Phase | PRP-1 |
| Format | Guard Cartridge Kit |
| Hardware Inner Diameter | 4.6 mm |
| Hardware Length | 20 mm |
| Particle Size | 12-20 µm |
| Packing Material Type | PSDVB |
| Guard Column Size | Semiprep/Preparative |
| Hardware Material | Stainless Steel |
| Max Pressure | 5,000 psi |
| Mobile Phase Limits | pH 1-13. 0-100% aqueous, organic modifier |
| Compound Classification | Organics |
| UOM | 1 holder, 1 cartridge |

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79121>