

PRP-1 Analytical Guard Cartridge Starter Kit (1 holder, 2 cartridges), Stainless Steel

Specification Sheet (Part/REF # 79447)



Organic compounds: small molecules (< 2,000 mw), oligonucleotides, pharmaceuticals, steroids, nucleic acids, vitamins, herbicides.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-1
Format	Guard Cartridge Kit
Hardware Inner Diameter	2.1 mm
Hardware Length	20 mm
Particle Size	10 µm
Packing Material Type	PSDVB
Guard Column Size	Analytical
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Organics
UOM	1 holder, 2 cartridges

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79447>

Spec. Version B