

# PRP-3 Analytical Guard Column Replacement Cartridges (5/pk), PEEK

Specification Sheet (Part/REF # 79395)



Protect your PRP-3 column.

## Product Specifications

<b>Separation Mode</b>	Reversed-Phase
<b>Stationary Phase</b>	PRP-3
<b>Format</b>	Replacement Guard Cartridge
<b>Hardware Inner Diameter</b>	3.0 mm
<b>Hardware Length</b>	8 mm
<b>Particle Size</b>	10 µm
<b>Packing Material Type</b>	PSDVB
<b>Guard Column Size</b>	Analytical
<b>Hardware Material</b>	PEEK
<b>Max Pressure</b>	5,000 psi
<b>Mobile Phase Limits</b>	pH 1-13. 0-100% aqueous, organic modifier
<b>Compound Classification</b>	Peptides/Proteins
<b>Frit Material</b>	Titanium
<b>UOM</b>	5/pk

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79395>

**Spec. Version B**