

PRP-3 Analytical Guard Column Starter Kit (1 holder, 2 cartridges), Stainless Steel

Specification Sheet (Part/REF # 79461)



Protect your PRP-3 column.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-3
Format	Guard Cartridge Kit
Hardware Inner Diameter	2.1 mm
Hardware Length	20 mm
Particle Size	10 µm
Packing Material Type	PSDVB
Guard Column Size	Analytical
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Peptides/Proteins
UOM	1 holder, 2 cartridges

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79461>

Spec. Version B