

PRP-h5 Semiprep/Prep Guard Column Starter Kit (1 holder, 1 cartridge), Stainless Steel

Specification Sheet (Part/REF # 79277)



Protect your PRP-h5 column.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-h5
Format	Guard Cartridge Kit
Hardware Inner Diameter	4.6 mm
Hardware Length	20 mm
Particle Size	12-20 μm
Packing Material Type	PSDVB
Particle Surface Area (m ² /cm ³)	277
Pore Volume (cm ³ /g)	0.526
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Peptides/Proteins
UOM	1 holder, 1 cartridge

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79277>

Spec. Version B