

PRP-C18

Semiprep/Preparative Guard Column Starter Kit (1 holder, 1 cartridge), Stainless Steel

Specification Sheet (Part/REF # 79689)



Protect your PRP-C18 column.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-C18
Format	Guard Cartridge Kit
Hardware Inner Diameter	4.6 mm
Hardware Length	20 mm
Particle Size	12-20 μ m
Packing Material Type	PS-DVB functionalized with C18
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Organics
UOM	1 holder, 1 cartridge

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79689>