

25 gauge, Large Hub RN Needle,
custom length (0.5 to 12 in), point
style 5 or AS, 6/PK

Specification Sheet (Part/REF # 7732-05)



Large Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 250 µL and larger, or SampleLock syringes. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

Product Specifications

Gauge	25 gauge
Point Style	5 or AS
Needle Length	Custom Needle Length
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	Large RN

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/needles/7732-05>

Spec. Version B