

25s gauge, Large Hub RN Needle,  
custom length (0.4 to 12 in), point  
style 2, 3, or 4, 6/PK

Specification Sheet (Part/REF # 7806-09)



Large Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 250 µL and larger, or SampleLock syringes. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

Product Specifications

Gauge	25s gauge
Point Style	2, 3, or 4
Needle Length	Custom Needle Length
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	Large RN

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/needles/7806-09>

Spec. Version B