

# 26 gauge, Large Hub RN Needle, 2 in, point style 3, 6/PK

Specification Sheet (Part/REF # 7780-02)



Large Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 250 µL and larger, or SampleLock syringes. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

## Product Specifications

Gauge	26 gauge
Point Style	3: Blunt
Needle Length	2.00 in (50.8 mm)
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	Large RN

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/needles/7780-02>

Spec. Version B