

**PEEK 0.062 inch OD/ 0.030 inch ID,
Large Hub RN Needle, 2 in, point
style 3, 6/PK**

Specification Sheet (Part/REF # 8650-01)



Kel-F and Removable Needles using standard HPLC grade PEEK tubing.

Product Specifications

Gauge	0.062 inch OD/ 0.030 inch ID
Point Style	3: Blunt
Needle Length	2.00 in (50.8 mm)
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	Large RN

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/needles/8650-01>

Spec. Version B