


# 0.2 mm Gel Loading, Small Hub RN Needle, 2.75 in, point style 3, 4/PK

Specification Sheet (Part/REF # 78633)



These are replacement needles for the adjustable width Gel Loading Syringe. The needles are available in 0.2 and 0.3 mm outer diameters. The needles are jacketed with a thicker gauge needle so only a short distance of the fine gauge needle is exposed.

## Product Specifications

<b>Gauge</b>	0.2 mm
<b>Point Style</b>	3
<b>Needle Length</b>	2.75 in (69.9 mm)
<b>Autoclavable</b>	Yes
<b>Gas Sterilizable</b>	Yes
<b>Sales Quantity</b>	4 needles per pack
<b>Hub</b>	Small RN

Specifications are subject to change without notice

**Spec. Version B**

<https://www.hamiltoncompany.com/laboratory-products/needles/78633>