


0.2 mm Gel Loading, Small Hub RN Needle, 2.75 in, point style 3, 4/PK

Specification Sheet (Part/REF # 78633)



These are replacement needles for the adjustable width Gel Loading Syringe. The needles are available in 0.2 and 0.3 mm outer diameters. The needles are jacketed with a thicker gauge needle so only a short distance of the fine gauge needle is exposed.

Product Specifications

Gauge	0.2 mm
Point Style	3: Blunt
Needle Length	2.75 in (69.9 mm)
Autoclavable	Yes
Gas Sterilizable	Yes
Sales Quantity	4 needles per pack
Hub	Small RN

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/needles/78633>

Spec. Version B