

# 0.3 mm Gel Loading, Small Hub RN Needle, 1 in, point style 3, 6/PK

Specification Sheet (Part/REF # 8651-02)



These are replacement needles for the fixed width Gel Loading Syringe. The needles are available in 0.2, 0.3, and 0.4 mm outer diameters. The 0.2 and 0.3 mm needles are jacketed with a thicker gauge needle so only a short distance of the fine gauge needle is exposed.

## Product Specifications

<b>Gauge</b>	0.3 mm
<b>Point Style</b>	3
<b>Needle Length</b>	1.00 in (25.4 mm)
<b>Autoclavable</b>	Yes
<b>Sales Quantity</b>	6 needles per pack
<b>Hub</b>	Small RN

Specifications are subject to change without notice

**Spec. Version B**

<https://www.hamiltoncompany.com/laboratory-products/needles/8651-02>