

# 22 gauge, Small Hub RN Needle, custom length (1.00 to 6.5 in), point style 3T, 6/PK

Specification Sheet (Part/REF # 8646-05)



Small Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 100  $\mu$ L and smaller. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

## Product Specifications

<b>Gauge</b>	22 gauge
<b>Point Style</b>	3T
<b>Needle Length</b>	Custom Needle Length
<b>Autoclavable</b>	Yes
<b>Sales Quantity</b>	6 needles per pack
<b>Hub</b>	Small RN

Specifications are subject to change without notice

**Spec. Version B**

<https://www.hamiltoncompany.com/laboratory-products/needles/8646-05>