

# 26s gauge, Small Hub RN Needle, custom length (1.00 to 6.5 in), point style 3T, 6/PK

Specification Sheet (Part/REF # 8646-02)



Small Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 100 µL and smaller. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

## Product Specifications

Gauge	26s gauge
Point Style	3T: Blunt PTFE Coated
Needle Length	Custom Needle Length
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	Small RN

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/needles/8646-02>

Spec. Version B