

26s gauge, Small Hub RN Needle, custom length (1.00 to 6.5 in), point style 3T, 6/PK

Specification Sheet (Part/REF # 8646-02)



Small Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 100 μ L and smaller. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

Product Specifications

Gauge	26s gauge
Point Style	3T
Needle Length	Custom Needle Length
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	Small RN

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/needles/8646-02>