

Adjustable Volume SoftGrip Pipette

10-100 µL

Specification Sheet (Part/REF # 55019-36)



Available in 0.2µL - 1 mL, SoftGrip’s isolated volume adjustment increases accuracy by removing parallax errors and preventing accidental volume adjustment, thanks to our next-generation dual indicator. Hamilton pipettors boast a slim design that allows for reaching the bottom of narrow tubes without the need for extended-length tips.

Product Specifications

Volume Range	10 to 100 µL
Channels	1
Color	Violet
Accuracy at 100 Percent Volume	0.60 ±%
Accuracy at 10% Volume	1.80 ±%
Precision at 100 Percent Volume	0.15 ±%
Precision at 10% Volume	0.70 ±%
Increments	0.2 µL
Autoclavable	Yes
Field Calibration	Yes

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/pipettes/55019-36>

Spec. Version B