

Adjustable Volume SoftGrip Pipette

2.5-25 µL

Specification Sheet (Part/REF # 55019-34)



Available in 0.2µL - 1 mL, SoftGrip’s isolated volume adjustment increases accuracy by removing parallax errors and preventing accidental volume adjustment, thanks to our next-generation dual indicator. Hamilton pipettors boast a slim design that allows for reaching the bottom of narrow tubes without the need for extended-length tips.

Product Specifications

Volume Range	2.5 to 25 µL
Channels	1
Color	Forest Green
Accuracy at 100 Percent Volume	0.80 ±%
Accuracy at 10% Volume	4.50 ±%
Precision at 100 Percent Volume	0.20 ±%
Precision at 10% Volume	1.50 ±%
Increments	0.02 µL
Autoclavable	Yes
Field Calibration	Yes

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/pipettes/55019-34>

Spec. Version B