

Adjustable Volume SoftGrip Pipette

30-300 µL

Specification Sheet (Part/REF # 55019-38)



Available in 0.2µL - 1 mL, SoftGrip’s isolated volume adjustment increases accuracy by removing parallax errors and preventing accidental volume adjustment, thanks to our next-generation dual indicator. Hamilton pipettors boast a slim design that allows for reaching the bottom of narrow tubes without the need for extended-length tips.

Product Specifications

Volume Range	30 to 300 µL
Channels	1
Color	Brick Red
Accuracy at 100 Percent Volume	0.40 ±%
Accuracy at 10% Volume	1.20 ±%
Precision at 100 Percent Volume	0.15 ±%
Precision at 10% Volume	0.40 ±%
Increments	0.2 µL
Autoclavable	Yes
Field Calibration	Yes

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/pipettes/55019-38>

Spec. Version B