

Fixed Volume SoftGrip Pipette 250 μ L

Specification Sheet (Part/REF # 55019-13)



SoftGrip fixed volume pipettes come with a slim ejection sleeve made of chemically resistant, high-strength polymer, as well as a universal calibration key to reduce calibration time. Hamilton fixed volume pipettes are available in 5 μ L - 1 mL.

Product Specifications

Volume Range	250 μ L Fixed
Channels	1
Color	Burnt Orange
Accuracy at 100 Percent Volume	0.40 \pm %
Accuracy at 10% Volume	N/A
Precision at 100 Percent Volume	0.18 \pm %
Precision at 10% Volume	N/A
Increments	N/A
Autoclavable	Yes
Field Calibration	Yes
Purity Level	Free of RNase/DNase, ATP, endotoxin, cytotoxic effects, and trace heavy metals

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/pipettes/55019-13>

Spec. Version B