

Fixed Volume SoftGrip Pipette

250 μ L

Specification Sheet (Part/REF # 55019-13)



The SoftGrip fixed volume pipettes are available in 5 μ L - 1 mL. The slim ejection sleeve is made from a high strength, chemically resistant polymer. The slim design allows the pipette to reach to the bottom of narrow tubes without the need for extended length tips.

Product Specifications

Volume Range	250 μ L Fixed
Channels	1
Color	Burnt Orange
Accuracy at 100 Percent Volume	0.40 \pm %
Accuracy at 10% Volume	N/A
Precision at 100 Percent Volume	0.18 \pm %
Precision at 10% Volume	N/A
Increments	N/A
Autoclavable	Yes
Field Calibration	Yes
Purity Level	Free of RNase/DNase, ATP, endotoxin, cytotoxic effects, and trace heavy metals

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/pipettes/55019-13>