

Fixed Volume SoftGrip Pipette

5 μ L

Specification Sheet (Part/REF # 55019-01)



The SoftGrip fixed volume pipettes are available in 5 μ L - 1 mL. The slim ejection sleeve is made from a high strength, chemically resistant polymer. The slim design allows the pipette to reach to the bottom of narrow tubes without the need for extended length tips.

Product Specifications

Volume Range	5 μ L Fixed
Channels	1
Color	Steel Blue
Accuracy at 100 Percent Volume	1.40 \pm %
Accuracy at 10% Volume	N/A
Precision at 100 Percent Volume	0.75 \pm %
Precision at 10% Volume	N/A
Increments	N/A
Autoclavable	Yes
Field Calibration	Yes

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/pipettes/55019-01>