

# Fixed Volume SoftGrip Pipette

## 50 $\mu$ L

Specification Sheet (Part/REF # 55019-07)



The SoftGrip fixed volume pipettes are available in 5  $\mu$ L - 1 mL. The slim ejection sleeve is made from a high strength, chemically resistant polymer. The slim design allows the pipette to reach to the bottom of narrow tubes without the need for extended length tips.

## Product Specifications

<b>Volume Range</b>	50 $\mu$ L Fixed
<b>Channels</b>	1
<b>Color</b>	Sandstone
<b>Accuracy at 100 Percent Volume</b>	0.60 $\pm$ %
<b>Accuracy at 10% Volume</b>	N/A
<b>Precision at 100 Percent Volume</b>	0.30 $\pm$ %
<b>Precision at 10% Volume</b>	N/A
<b>Increments</b>	N/A
<b>Autoclavable</b>	Yes
<b>Field Calibration</b>	Yes

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/pipettes/55019-07>

**Spec. Version B**