

# Fixed Volume SoftGrip Pipette 50 µL

Specification Sheet (Part/REF # 55019-07)



SoftGrip fixed volume pipettes come with a slim ejection sleeve made of chemically resistant, high-strength polymer, as well as a universal calibration key to reduce calibration time. Hamilton fixed volume pipettes are available in 5 µL - 1 mL.

## Product Specifications

Volume Range	50 µL Fixed
Channels	1
Color	Sandstone
Accuracy at 100 Percent Volume	0.60 ±%
Accuracy at 10% Volume	N/A
Precision at 100 Percent Volume	0.30 ±%
Precision at 10% Volume	N/A
Increments	N/A
Autoclavable	Yes
Field Calibration	Yes

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/pipettes/55019-07>

Spec. Version B