

# Fixed Volume SoftGrip Pipette 500 $\mu$ L

Specification Sheet (Part/REF # 55019-17)



SoftGrip fixed volume pipettes come with a slim ejection sleeve made of chemically resistant, high-strength polymer, as well as a universal calibration key to reduce calibration time. Hamilton fixed volume pipettes are available in 5  $\mu$ L - 1 mL.

## Product Specifications

Volume Range	500 $\mu$ L Fixed
Channels	1
Color	Olive Green
Accuracy at 100 Percent Volume	0.40 $\pm$ %
Accuracy at 10% Volume	N/A
Precision at 100 Percent Volume	0.18 $\pm$ %
Precision at 10% Volume	N/A
Increments	N/A
Autoclavable	Yes
Field Calibration	Yes

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/pipettes/55019-17>

Spec. Version B