

# Fixed Volume SoftGrip Pipette

## 500 µL

Specification Sheet (Part/REF # 55019-17)



The SoftGrip fixed volume pipettes are available in 5 µL - 1 mL. The slim ejection sleeve is made from a high strength, chemically resistant polymer. The slim design allows the pipette to reach to the bottom of narrow tubes without the need for extended length tips.

## Product Specifications

<b>Volume Range</b>	500 µL Fixed
<b>Channels</b>	1
<b>Color</b>	Olive Green
<b>Accuracy at 100 Percent Volume</b>	0.40 ±%
<b>Accuracy at 10% Volume</b>	N/A
<b>Precision at 100 Percent Volume</b>	0.18 ±%
<b>Precision at 10% Volume</b>	N/A
<b>Increments</b>	N/A
<b>Autoclavable</b>	Yes
<b>Field Calibration</b>	Yes

Specifications are subject to change without notice

**Spec. Version B**

<https://www.hamiltoncompany.com/laboratory-products/pipettes/55019-17>