

12-Channel Adjustable Volume SoftGrip Pipette 5-50 µL

Specification Sheet (Part/REF # 55019-54)



Hamilton multi-channel pipettes come with either 8 or 12 channels and are available for both 5 µL - 50 µL or 30 µL - 300 µL volumes. The movable lower body and unique long-tip ejector allow for full 360° rotation to provide maximum comfort for both right-handed and left-handed users.

Product Specifications

Volume Range	5 µL to 50 µL
Channels	12
Color	Sandstone
Accuracy at 100 Percent Volume	1.00 ±%
Accuracy at 10% Volume	5.00 ±%
Precision at 100 Percent Volume	0.50 ±%
Precision at 10% Volume	2.00 ±%
Increments	0.02 µL
Autoclavable	Yes
Field Calibration	Yes

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/pipettes/55019-54>

Spec. Version B