


# 8-Channel Adjustable Volume SoftGrip Pipette 30-300 µL

Specification Sheet (Part/REF # 55019-52)



Hamilton multi-channel pipettes come with either 8 or 12 channels and are available for both 5 µL - 50 µL or 30 µL - 300 µL volumes. The movable lower body and unique long-tip ejector allow for full 360° rotation to provide maximum comfort for both right-handed and left-handed users.

## Product Specifications

Volume Range	30 µL to 300 µL
Channels	8
Color	Brick Red
Accuracy at 100 Percent Volume	0.50 ±%
Accuracy at 10% Volume	1.20 ±%
Precision at 100 Percent Volume	0.30 ±%
Precision at 10% Volume	0.70 ±%
Increments	0.2 µL
Autoclavable	Yes
Field Calibration	Yes

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/pipettes/55019-52>

Spec. Version B