

8-Channel Adjustable Volume SoftGrip Pipette 5-50 µL

Specification Sheet (Part/REF # 55019-50)



Hamilton multi-channel pipettes come with either 8 or 12 channels and are available for both 5 µL - 50 µL or 30 µL - 300 µL volumes. The movable lower body and unique long-tip ejector allow for full 360° rotation to provide maximum comfort for both right-handed and left-handed users.

Product Specifications

| | |
|---------------------------------|---------------|
| Volume Range | 5 µL to 50 µL |
| Channels | 8 |
| Color | Sandstone |
| Accuracy at 100 Percent Volume | 0.80 ±% |
| Accuracy at 10% Volume | 2.00 ±% |
| Precision at 100 Percent Volume | 0.35 ±% |
| Precision at 10% Volume | 2.00 ±% |
| Increments | 0.02 µL |
| Autoclavable | Yes |
| Field Calibration | Yes |

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/pipettes/55019-50>

Spec. Version B