


8-Channel Adjustable Volume SoftGrip Pipette 5-50 μ L

Specification Sheet (Part/REF # 55019-50)



Hamilton multi-channel pipettes come with either 8 or 12 channels and are available for both 5 μ L - 50 μ L or 30 μ L - 300 μ L volumes. The movable lower body and unique long-tip ejector allow for full 360° rotation to provide maximum comfort for both right-handed and left-handed users.

Product Specifications

Volume Range	5 μ L to 50 μ L
Channels	8
Color	Sandstone
Accuracy at 100 Percent Volume	0.80 \pm %
Accuracy at 10% Volume	2.00 \pm %
Precision at 100 Percent Volume	0.35 \pm %
Precision at 10% Volume	2.00 \pm %
Increments	0.02 μ L
Autoclavable	Yes
Field Calibration	Yes

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/pipettes/55019-50>

Spec. Version B