

# SoftGrip Multi-Channel Stand

Specification Sheet (Part/REF # 34195)



Hamilton's lab pipettes and pipettors work with a variety of accessories, including carousels, SoftGrip stands, and calibration equipment, designed to help you maximize quality, comfort, and precision.

## Product Specifications

Pipette Type	Manual (Multi-Channels only)
Specifications are subject to change without notice	
<a href="https://www.hamiltoncompany.com/laboratory-products/pipettes/34195">https://www.hamiltoncompany.com/laboratory-products/pipettes/34195</a>	
Spec. Version B	