

1000 μ L, Bulk, Low Retention, 1000 tips/bag

Specification Sheet (Part/REF # 68010-20)

Hamilton delivers performance, value, and choice with two high quality tip families, plus multiple formats and packaging options. Formats include standard, non-filter, extended length, hydrophobic filter barrier, non-sterile, pre-sterilized, and macrovolume.

Product Specifications

Disposable Tip Volume	1000 μ L
Tip Length	Standard
Packaging	Bulk
Quantity Per Each	1000
Filtered	No
Low Retention	Yes
Pre-Sterilized	No

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/pipettes/68010-20>