

1000 μ L, Racked, Filter, Pre-Sterilized, 96 tips/rack

Specification Sheet (Part/REF # 68010-13)

Hamilton delivers performance, value, and choice with two high quality tip families, plus multiple formats and packaging options. Formats include standard, non-filter, extended length, hydrophobic filter barrier, non-sterile, pre-sterilized, and macrovolume.

Product Specifications

Disposable Tip Volume	1000 μ L
Tip Length	Standard
Packaging	Racked
Quantity Per Each	576
Filtered	Yes
Low Retention	No
Pre-Sterilized	Yes

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/pipettes/68010-13>