

# 10 µL, Bulk, Low Retention, 1000 tips/bag

Specification Sheet (Part/REF # 68007-20)



Hamilton pipette tip families deliver value, choice, and performance across a variety of formats and packaging options. Formats include standard, non-filter, extended length, hydrophobic filter barrier, non-sterile, sterilized, and macro volume.

## Product Specifications

Disposable Tip Volume	10 µL
Tip Length	Standard
Packaging	Bulk
Quantity Per Each	1000
Filtered	No
Low Retention	Yes
Sterile	No
Purity Level	Free of RNase/DNase, ATP, endotoxin, cytotoxic effects, and trace heavy metals
Recommended Pipette Volume Ranges	2 µL, 5 µL, and 10 µL Pipettes

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/pipettes/68007-20>

Spec. Version B