

# 10 µL, Racked, Low Retention, 96 tips/rack

Specification Sheet (Part/REF # 68007-21)



Hamilton pipette tip families deliver value, choice, and performance across a variety of formats and packaging options. Formats include standard, non-filter, extended length, hydrophobic filter barrier, non-sterile, sterilized, and macro volume.

## Product Specifications

Disposable Tip Volume	10 µL
Tip Length	Standard
Packaging	Racked
Quantity Per Each	1152
Filtered	No
Low Retention	Yes
Sterile	No
Purity Level	Free of RNase/DNase, ATP, endotoxin, cytotoxic effects, and trace heavy metals
Recomended Pipette Volume Ranges	2 µL, 5 µL, and 10 µL Pipettes

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/pipettes/68007-21>

Spec. Version B