

# 10 $\mu$ L, Racked, Low Retention, 96 tips/rack

Specification Sheet (Part/REF # 68007-21)

Hamilton delivers performance, value, and choice with two high quality tip families, plus multiple formats and packaging options. Formats include standard, non-filter, extended length, hydrophobic filter barrier, non-sterile, pre-sterilized, and macrovolume.

## Product Specifications

|                              |  |
|------------------------------|--|
| <b>Disposable Tip Volume</b> | 10 $\mu$ L   |
| <b>Tip Length</b>            | Standard   |
| <b>Packaging</b>             | Racked   |
| <b>Quantity Per Each</b>     | 1152   |
| <b>Filtered</b>              | No   |
| <b>Low Retention</b>         | Yes  |
| <b>Pre-Sterilized</b>        | No   |
| <b>Purity Level</b>          | Free of RNase/DNase, ATP, endotoxin, cytotoxic effects, and trace heavy metals |

Specifications are subject to change without notice

**Spec. Version B**

<https://www.hamiltoncompany.com/laboratory-products/pipettes/68007-21>