

10 μ L, Racked, Extended Length, Filter, Pre-Sterilized, Low Retention, 96 tips/rack

Specification Sheet (Part/REF # 68007-26)

Hamilton delivers performance, value, and choice with two high quality tip families, plus multiple formats and packaging options. Formats include standard, non-filter, extended length, hydrophobic filter barrier, non-sterile, pre-sterilized, and macrovolume.

Product Specifications

Disposable Tip Volume	10 μ L
Tip Length	Extended
Packaging	Racked
Quantity Per Each	960
Filtered	Yes
Low Retention	Yes
Pre-Sterilized	Yes
Purity Level	Free of RNase/DNase, ATP, endotoxin, cytotoxic effects, and trace heavy metals

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/pipettes/68007-26>