

10 µL, Racked, Extended Length, Filter, Sterile, Low Retention, 96 tips/rack

Specification Sheet (Part/REF # 68007-26)



Hamilton pipette tip families deliver value, choice, and performance across a variety of formats and packaging options. Formats include standard, non-filter, extended length, hydrophobic filter barrier, non-sterile, sterilized, and macro volume.

Product Specifications

| | |
|----------------------------------|--|
| Disposable Tip Volume | 10 µL |
| Tip Length | Extended |
| Packaging | Racked |
| Quantity Per Each | 960 |
| Filtered | Yes |
| Low Retention | Yes |
| Sterile | Yes |
| Purity Level | Free of RNase/DNase, ATP, endotoxin, cytotoxic effects, and trace heavy metals |
| Recomended Pipette Volume Ranges | 2 µL, 5 µL, and 10 µL Pipettes |

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/pipettes/68007-26>

Spec. Version B