

# 200 µL, Racked, Filter, Low Retention, Sterile, 96 tips/rack

Specification Sheet (Part/REF # 68008-23)



Hamilton pipette tip families deliver value, choice, and performance across a variety of formats and packaging options. Formats include standard, non-filter, extended length, hydrophobic filter barrier, non-sterile, sterilized, and macro volume.

## Product Specifications

Disposable Tip Volume	200 µL
Tip Length	Standard
Packaging	Racked
Quantity Per Each	960
Filtered	Yes
Low Retention	Yes
Sterile	Yes
Purity Level	Free of RNase/DNase, ATP, endotoxin, cytotoxic effects, and trace heavy metals
Recomended Pipette Volume Ranges	20 µL, 25 µL, 50 µL, 100 µL, and 200 µL Pipettes

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/pipettes/68008-23>

Spec. Version B