

# 200 $\mu$ L, Racked, Low Retention, Pre-Sterilized, 96 tips/rack

Specification Sheet (Part/REF # 68008-22)

Hamilton delivers performance, value, and choice with two high quality tip families, plus multiple formats and packaging options. Formats include standard, non-filter, extended length, hydrophobic filter barrier, non-sterile, pre-sterilized, and macrovolume.

## Product Specifications

<b>Disposable Tip Volume</b>	200 $\mu$ L
<b>Tip Length</b>	Standard
<b>Packaging</b>	Racked
<b>Quantity Per Each</b>	960
<b>Filtered</b>	No
<b>Low Retention</b>	Yes
<b>Pre-Sterilized</b>	Yes
<b>Purity Level</b>	Free of RNase/DNase, ATP, endotoxin, cytotoxic effects, and trace heavy metals

Specifications are subject to change without notice

**Spec. Version B**

<https://www.hamiltoncompany.com/laboratory-products/pipettes/68008-22>